



# BP+ User's Manual



Central Blood Pressure Measuring Device



© Copyright Uscom Limited 2009-2014. All rights reserved.

No part of this document may be reproduced or transmitted in any form or by any means, electronic, mechanical, photocopying, recording, or otherwise, without the prior written permission of Uscom Limited. Product specifications may change without prior notice.

CardioScope, Pulsecor and BP+ are trademarks of Uscom Limited.

#### Disclaimers and Limitations of Liability

The information contained in this document, including but not limited to any product specifications, is subject to change without notice. Uscom makes no representations or warranties of any kind concerning the contents of this document. Under no circumstances will Uscom be held liable for any loss or other damages pertaining to the use of this document.

 <p>Uscom Limited Level 7, 10 Loftus St Sydney 2000 NSW Australia  T: +61 2 92474144 F: +61 2 92478157  E: <a href="mailto:info@uscom.com.au">info@uscom.com.au</a> W: <a href="http://www.uscom.com.au">www.uscom.com.au</a></p>	 <p>Emergo Eurpoe Molenstraat 15 2513 BH, The Hague The Netherlands  T: +31 70 345 8570 F: +31 70 346 7299  W: <a href="http://www.emergogroup.com">www.emergogroup.com</a></p>
--	--

#### BP+ User Manual (R7)

Part Number: 0760024-00

Pertains To: BP+ (BPPLUS-R7)

Version: 5.0, 2014-02



# Table of Contents

How to Use this Manual.....	4
Glossary .....	5
Symbols .....	6
Safety Information .....	8
1 Product Overview .....	10
1.1 Intended Use .....	10
1.2 Setting Up .....	10
1.3 Supplying Power .....	10
1.3.1 Turning Device On .....	10
1.3.2 Turning Device Off .....	11
1.3.3 Power On Self-Tests .....	11
1.4 Connecting a Cuff to BP+.....	12
1.5 Saving Measurements .....	12
2 Common Procedures.....	13
2.1 Preparing for a measurement .....	13
2.2 Placing the cuff .....	13
2.3 Taking a measurement .....	13
2.4 Displaying the Result.....	14
2.5 Cancelling the measurement .....	15
2.6 Interpreting Signal Quality .....	15
2.7 Displaying Pulse Wave Analysis.....	16
2.8 Settings Menu .....	16
2.9 Setting Inflation Pressure .....	17
2.10 Setting Date and Time .....	18
2.11 Selecting Language.....	18
2.12 Report Printing & Report Style .....	19
2.13 Reviewing saved measurements .....	20
3 Theory of Operation .....	22
3.1 Physiology of wave reflection .....	22
3.2 Blood Pressure Measurement .....	22
3.3 Display of Pulse Waveform and Rhythm strip .....	23
3.4 Central Blood Pressure .....	23
3.5 Pulse Wave Parameters Calculations .....	24
4 Maintenance and Troubleshooting .....	25
4.1 Servicing .....	25
4.2 Routine Maintenance .....	25
4.3 Cleaning.....	25
4.4 Disposal .....	25
4.5 Frequently Asked Questions .....	25
4.6 Error Message .....	27
4.6.1 Main Screen Message .....	27
4.6.2 Memory Mode Message.....	28
4.6.3 SD Card Icon Message .....	29
5 Specifications.....	30
5.1 Electromagnetic Compatibility .....	31
6 Accessories and Spare Parts .....	32

## How to Use this Manual

This manual uses some conventions intended to make it easier to understand.



Warning statements convey safety information to mitigate against the potential for death or injury and are displayed in a red box with a warning symbol to the left.



Caution statements convey information about the potential to damage equipment, produce inaccurate data or invalidate a procedure and are displayed in a dashed box with a caution symbol to the left.



Statements conveying helpful hints and notes are displayed with an information symbol to the left.

*Individual instructions* are displayed in blue italics with accompanying text.

Text displayed on the device screen is drawn enclosed in a box.














## Glossary

The following terms and abbreviations are used in this manual.

Term	Description
Sys BP	Systolic Blood Pressure
Dia BP	Diastolic Blood Pressure
Central Sys BP	Central Systolic Blood Pressure
Central Dia BP	Central Diastolic Blood Pressure
AI	Augmentation Index
PR	Pulse Rate
Incident wave	A wave propagating directly from the heart
NIBP	Non-invasive blood pressure
PC	Personal computer
POST	Power On Self-Test
Reflected wave	A wave propagating from a reflection site
SD	Secure Data, type of non-volatile removable storage card
Slot	A place to store a single measurement
Pulse wave	Pertaining to something above systolic, for example a signal recorded at a cuff pressure above systolic pressure





## Symbols

The following symbols appear on the BP+ device and in accompanying labelling and documentation.

Symbol	Description
	Class II equipment (IEC 60601-1)
	Type BF applied part (IEC 60601-1) - Blood Pressure Cuff
	Direct current
	Warning, consult safety information
	Consult instructions for use
	Information
	Caution
IOIOI	Serial interface
	USB interface
	Do not dispose of the BP+ as unsorted municipal waste
	Manufacturer
	Date of manufacture
	Device serial number
	Device model number





**Table 1: Symbols on BP+**

The following symbols when used in this manual refer to hardware buttons on the device.

Symbol	Button Name	Functions
	Start	Start measurement Confirm selection Cancel measurement in progress
	Next	Move to the next item Enter memory mode Cancel measurement in progress
	Return	Exit a menu Display additional measurement information Cancel measurement in progress
	Previous	Move to previous item Enter inflation target mode Cancel measurement in progress

**Table 2: Buttons on BP+**

The following symbols when used in this manual refer to icons that appear on the device LCD Screen.

Symbol	Description
	Adult mode
	Paediatric mode
	Secure Data (SD) Card inserted and ready
	Secure Data (SD) Card inserted but BP+ cannot use the card

**Table 3: BP+ screen icons**

## Safety Information



### **Electric Shock Hazards:**

- The BP+ may be isolated from the mains supply by disconnecting the power supply from the mains supply outlet or by disconnecting the power supply from the monitor.
- Do not open the enclosure of the BP+ or power supply unit, especially while they are connected to an AC supply outlet.
- Disconnect the BP+ and power supply unit from the AC supply outlet before cleaning. Do not use liquid or spray detergents.
- Avoid ingress of liquid into any part of the BP+ or power supply unit. Do not submerge any component in liquid. This may cause fire or electrical shock.

**Explosion Hazard:** Do not use the BP+ in a flammable atmosphere or where concentrations of flammable anaesthetics may occur.

### **Power Supply:**

- Ensure the AC supply voltage is correct and the correct mains power adaptor is securely fitted in place before connecting the BP+ and mains power supply to the AC supply outlet.
- Use only the mains power adaptor supplied with the BP+ or a genuine replacement part (GlobTek Inc, model GTM41060-1706).
- This product is intended for indoor use only.
- Protect from excessive force or shock.
- Do not pull output plug with excessive force.

**Accessories and Equipment:** Use of accessories or equipment not approved by Uscom, or not complying with safety standards equivalent to those met by the BP+, may lead to a reduced level of safety of the resulting system or failure of the BP+ to operate correctly.

**Connection to Personal Computer:** Equipment connected to the serial or USB port of the BP+ must be certified to applicable IEC standards (IEC 60950 for data processing equipment and IEC 60601-1 for medical electrical equipment). All configurations must comply with the medical electrical equipment standard IEC 60601-1 Clause 16. Anyone who connects equipment to the BP+ configures a medical system, and is responsible for ensuring that the system complies with the requirements of the medical electrical equipment standard IEC 60601-1 Clause 16.

**Device Modification:** No modification of the BP+ is allowed.





**Removal from Use:** If any of the following situations occur, stop using the BP+:

- The AC supply cable or attachment plug is damaged.
- The device has been exposed to excessive moisture.
- The device has been dropped and damaged or shows obvious signs of breakage.
- The device display remains blank when turned on.

**Contact with Patient:** Do not touch any exposed metal parts on the BP+ (in particular interface connectors) while simultaneously touching a patient.

Too frequent measurements can cause injury to the patient due to blood flow interference.

Measurement in an arm where intravascular access or therapy or an arterio-venous shunt is present could result in injury to the patient.

Pressure applied during measurement may cause discomfort or injury where the patient has had a mastectomy on the side the cuff is applied.

**Operator Attendance:** The operator must be in continual attendance while using the BP+.



**Caution:** Do not clean the BP+ with concentrated bleach, corrosive chemicals or abrasive cleaning compounds.

**Electromagnetic Compatibility:** The BP+ complies with the electromagnetic compatibility requirements of IEC 60601-1-2. Operation of the device may affect, or be affected by, nearby equipment due to the effects of electromagnetic interference. If this happens:

- Increase the separation between the BP+ and the other device.
- Connect the devices to AC supply outlets on separate circuit branches.
- Refer to section 5.1 for further compliance information and advice relating to electromagnetic interference.

The use of a mains power adapter not approved by Uscom may result in increased electromagnetic emissions or decreased electromagnetic immunity of the BP+.

# 1 Product Overview

Congratulations on your acquisition of the Uscom BP+. Designed for ease of use, accuracy and repeatability, the BP+ will easily integrate into many different clinical and research environments. The BP+ measures brachial and central blood pressure, and other cardiac and arterial parameters, using the cuff oscillometric method.



BP+ central blood pressure measurements have been validated on patients who are eighteen (18) and older.

## 1.1 Intended Use

BP+ is a non-invasive, compact standalone measurement device that automatically measures systolic and diastolic pressure, and pulse rate in adult and paediatric patients. BP+ also provides non-invasive central (aortic) systolic and diastolic blood pressure, augmentation index, pulse waveform and other parameters intended for use in adult patients.

BP+ performs measurements using a conventional oscillometric method via a brachial cuff on the upper arm.

The device is intended to be used under supervision by qualified healthcare personnel. The device is not intended for operating theatre, intensive care, or continuous patient monitoring use.



**Caution:** BP+ should not be used on neonatal subjects. It is not intended for measurement of blood pressure in subjects who are pregnant.

**Caution:** Federal (USA) law restricts this device to sale by or on order of a physician.

## 1.2 Setting Up

To start using the device you will need to supply power, connect an inflatable cuff, and insert an SD Card if you wish to save the measurement data.

## 1.3 Supplying Power

The BP+ requires DC power input to the device. The mains power adapter provided with the monitor can provide the required DC power from most international mains power outlets from 100 to 240 VAC and 50 to 60 Hz mains input.

The mains power adapter is supplied with interchangeable blades. *Attach the blades that match your wall sockets.*

Insert the plug at the end of the cable attached to the AC wall adapter into the device in the socket indicated as the External DC input socket labelled **===**.



**Use only the mains power adapter supplied with the BP+ or a genuine replacement part (GlobTek Inc, model GTM41060-1706).**

### 1.3.1 Turning Device On

The device will turn on automatically when the AC power is attached.

### 1.3.2 Turning Device Off

To turn off the device, remove the AC power.



Turning the device off will release pressure in the cuff. A more effective way to quickly release pressure from a subject's arm is to remove the cuff from their arm.

### 1.3.3 Power On Self-Tests

Once the device is turned on, the device will load the BP+ application into the device internal memory. During this time it will display the following screen.




Figure 1: Boot up screen

BP+ will begin execution once loaded. The BP+ initial screen will report the version of BP+ and start a power on self-test (POST). An example initial screen is shown in Figure 2.



Figure 2: Start-up POST Screen



The version of software in the device may differ from what is displayed here.

Once the POST is complete, the device will clear the display and present the Device Ready screen as shown in Figure 3.




Figure 3: Ready Screen

### 1.4 Connecting a Cuff to BP+

In order to take a measurement, you will need to connect an inflatable blood pressure cuff using the extension hose supplied with the device.

*To connect the extension hose* push the air fittings together until you hear and feel a positive “click”.

*To disconnect the extension hose*, uncouple the fittings by fully depressing the thumb button and pulling the fittings apart gently.

*To connect a cuff*, push the extension hose fitting into the circular orifice in the cuff. You should hear and feel a positive “click”.

*To disconnect a cuff*, squeeze the sides of the extension hose fitting and remove it gently from the cuff.



The BP+ has been developed using specific cuff designs. Using alternate cuff designs may significantly alter the accuracy of measurement.



Avoid compression or kinking the air hoses connecting the device to the cuff. Such actions, particularly during a measurement, may result in prolonged inflation of the cuff which cannot be resolved by the device’s built-in safeguards. Prolonged over-inflation of a cuff may result in harm to the patient.

The operator should ensure that application of the cuff does not result in prolonged impairment of the circulation of blood.

Do not connect cuffs with luer lock connectors to Intravenous System (IV).


### 1.5 Saving Measurements



The BP+ can store measurements if an SD Card is inserted. The BP+ supports full size SD Cards with a storage capacity up to 2GB.

Insert the SD Card with the logo facing up to the user.

Note: Measurements will not be stored if a compatible SD Card is not inserted into the BP+. Look for the SD Card icon on the top left of the screen. If the icon is not present, then no SD Card is inserted. If the icon has a red cross through it, then the BP+ does not recognise the SD Card.



## 2 Common Procedures

The BP+ is designed for simple operation through the four buttons on the front of the device. These buttons are described in Table 2.

### 2.1 Preparing for a measurement

The entire measurement procedure usually takes less than 90 seconds. However, it is recommended that preparations for taking a measurement follow similar good practice procedures to those for taking non-invasive blood pressure.

*Ensure the subject is sitting* comfortably upright with adequate back support. The subject's legs should be uncrossed and feet should be flat on the floor.

*Instruct the subject* to relax, remain still, breathe normally, cease talking and avoid any movement during measurement. It is recommended that the subject is at rest for at least five minutes before the first measurement and to repeat the measurement at least three times, recording the median of the measurements.

### 2.2 Placing the cuff

*Apply the deflated cuff* to the subject's left upper arm. Choose a cuff of the correct size, as indicated by the range markings on the cuff. The cuff should be wrapped around the arm firmly, but not tightly. The artery marker on the cuff where the hose enters the cuff should be placed over the brachial artery. Avoid placing the cuff over thick clothing. Avoid bunching up clothing above the cuff, which may partially occlude the artery.



To reduce the risk of cross-infection and further abrasion, do not apply the cuff to broken skin.

*Support the subject's forearm* on a stable, horizontal surface such as a table top, at a height such that the cuff is at the approximate level of the subject's heart and held away from the subject's torso. Generally, a palm-down position is found to be more comfortable.

### 2.3 Taking a measurement

Press Start  button while the main screen is displayed as in Figure 4:





Figure 4: Main measurement screen

There may be a short pause while the device performs a pre-measurement check. The device will then begin inflating the cuff. During this time the status area will display **Measuring Blood Pressure** and the current cuff pressure will be displayed.




**Figure 5: Measuring screen**

If your model of BP+ supports “Measurement on Inflation”, the device will attempt to determine the systolic, mean and diastolic pressures and pulse rate as the cuff is inflated. If the model does not support measurement on inflation, or the device is unable to determine the blood pressure during inflation, once the pre-set target pressure is reached, the device will gradually deflate the cuff at a controlled rate of approximately 3 mmHg/s. During this time, the device will attempt to determine the systolic, mean and diastolic pressures, as well as pulse rate as the cuff is deflated. If the initial inflation pressure is less than the systolic pressure, the device will re-inflate the cuff to a higher pressure before the controlled deflation.

If your model of BP+ supports “Measurement on Inflation”, and it successfully determines the blood pressure during inflation, the device will directly inflate the cuff to the pulse wave pressure. Otherwise, once the blood pressure and pulse rate has been determined by deflating the cuff, the device will vent the remaining pressure in the cuff to a pressure below 5 mmHg, after which it will re-inflate the cuff to approximately 30 mmHg above the measured systolic pressure. This pulse wave pressure will be maintained for around 12 seconds. During this time, the status area will display Measuring Pulse Wave. After the pulse wave measurement is finished, the cuff pressure will be completely released. This marks the completion of the entire measurement cycle. The BP+ will then process the acquired data and display the result on the screen.

#### **2.4 Displaying the Result**

Once the measurement is completed, the monitor will display all results under the respective labels as follows:



**Figure 6: BP Result Screen**

10 seconds of the pulse wave signal are displayed along the bottom of the display. Vertical red lines indicate the start of each pulse. This display can provide a visual indication of any arrhythmias and inter-beat variability (such as that caused by respiratory load).



Measurement results may not be accurate in the presence of arrhythmias. Arrhythmias may be apparent on the rhythm strip.

The main screen will also display the average pulse wave pulse.

## 2.5 Cancelling the measurement

*Press any button* during a measurement to cancel the measurement and deflate the cuff.

Note: While the device is calculating results it will not respond to button presses.

## 2.6 Interpreting Signal Quality


The signal quality of the pulse wave measurements will be measured by the device and expressed as a Signal-to-Noise Ratio (SNR) on a logarithmic scale (dB). An increase in the SNR by 3 indicates a tenfold increase in the signal quality. The SNR acceptability is indicated on the results screen as a text status message and can vary from Invalid, Poor, Acceptable, Good, through to Excellent as shown in Table 4: Signal quality classification. The numeric value of SNR is displayed on the Pulse Wave Analysis screen as shown in Figure 7.

Measurements with a SNR below 6 (Poor or Invalid) should be discarded. The device will not display the augmentation index and central pressures in these circumstances. It is recommended to repeat the measurement and to confirm the patient is remaining still and not talking during the measurement. (Patients often move when talking.)

Signal-to-Noise Ratio (dB)	Signal Quality Classification	Colour on Results screen
SNR < 0	Invalid	Red
0 ≤ SNR < 6	Poor	Red
6 ≤ SNR < 9	Acceptable	Yellow
9 ≤ SNR < 12	Good	Green
12 ≤ SNR	Excellent	Green

**Table 4: Signal quality classification**

## 2.7 Displaying Pulse Wave Analysis

Press  button from the main screen. The monitor will display the Pulse Wave Analysis screen as follows:

Pulse Wave Analysis	
PRV: 0.0	PR (bpm): 80.0
SNR: 18.6	TR: 0.16
SEP: 0.26	RWTT: 0.14
PP: 5.1	DPP: 32.0 %
dP/dt(max): 98.89	Date: 2012/06/20
dP/dt(artery): 919.83	Time: 14:20:49

Figure 7: Pulse Wave Analysis Screen



## 2.8 Settings Menu

The BP+ provides a Settings menu for configuring different functions of the monitor. The setting options include: Inflation Pressure\*, Report Printing\*, Report Style\*, Select Language\* and Set Date and Time.

Settings
<b>1. Exit Settings Menu</b>
<b>2. Inflation Pressure *: Adult - 150 mmHg</b>
<b>3. Report Printing *: Disabled</b>
<b>4. Report Style *: Short</b>
<b>5. Select Language *: English</b>
<b>6. Set Date and Time</b>

\* Reset will be required if changed

Figure 8: Settings Menu Screen

Changing some setting options require a device reset as indicated by an asterisk (\*). Whenever a change is made to an asterisk (\*) setting option, the first option – Exit Settings Menu on the Settings menu will be changed to Save Settings and Reset. Therefore, you can either select Save Settings and Reset and then press  to reset the monitor, or you can press  to reset the monitor directly.






Figure 9: Settings Menu Screen – Save Settings and Reset

## 2.9 Setting Inflation Pressure

By default, the BP+ will initially inflate the cuff to 150 mmHg. If the actual systolic pressure is higher than this, then the device may automatically re-inflate additional times to measure the systolic pressure. If you have an estimate of the systolic pressure then you may wish to set the initial inflation target to avoid these additional inflation cycles. The monitor supports adult and paediatric mode inflation target ranges.

Press . The Settings menu will be displayed. Select the **Inflation Pressure** option and press . The Inflation Target menu will be displayed with the current inflation target highlighted. Press to confirm your selection or use the arrow keys to select another value. You will be returned to the Settings screen, and a reset is required - refer to **2.8 Settings Menu** for device reset. If you do not wish to confirm your selection, press .






Figure 10: Inflation Target Setup Screen


**Note** - If the subject's actual systolic pressure is lower than 150 mmHg you may reduce measurement time and discomfort by reducing the initial inflation target.

If a lower range of inflation targets would be more appropriate for your subject, you can change the list of available inflation targets by selecting the **Change to Paediatric Mode** option and pressing to confirm your selection. To return to a higher range of inflation targets, select **Change to Adult Mode** and press .

## 2.10 Setting Date and Time





The BP+ provides a **Set Date and Time** menu for setting the date and time values on the monitor.

Press . The Settings menu will be displayed. Select the **Set Date and Time** option and press . The date time setting menu will be displayed.





The screen is divided into two sections: 'Time' and 'Date'. The 'Time' section has three input fields: 'Hour' (13), 'Minute' (00), and 'Second' (31). The 'Date' section has three input fields: 'Year' (2012), 'Month' (01), and 'Day' (01). The fields are separated by colons and slashes.

Figure 11: Date Time Setting Screen

Press  to highlight a different date time element. Then use  and  to change the value of the highlighted date time element. The monitor will automatically store the changed value and does not require any reset of the monitor. To return to the main screen, press .

## 2.11 Selecting Language


By default, the language displayed on the monitor is English. The BP+ monitor may supports different language display for some countries.


Press . The Settings menu will be displayed. Then go to select the **Select Language** option and press . The Select Language menu will be displayed with the current language highlighted.



The screen is titled 'Select Language'. It has a list of language options: English, Deutsch, italiano, and español. The 'English' option is highlighted with a blue background.

Figure 12: Select Language Screen

Use the arrow keys to select the desired language. If you do not wish to confirm your selection, press  to return to the Settings menu screen. Otherwise press



 *to confirm your selection.* You will be returned to the Settings menu screen and a reset is required - refer to **2.8 Settings Menu** for device reset.



Some versions of monitor might only provide one language option for some countries. Please contact Uscom for more information using the contact information at the front of this manual.

## 2.12 Report Printing & Report Style

The BP+ monitor can be connected to a thermal printer to print the measurement report. By default, the report printing function is disabled. You can enable this function by changing the settings on the Settings Menu.

Press . The Settings menu will be displayed. Select the **Report Printing** option, *press  to toggle this option to Enabled.* Then reset the monitor as described in **2.8 Settings Menu**.





Figure 13: Enabling Report Printing Screen

A report guidance table can be printed along with the measurement report. To include the report guidance table to the report printing, select the **Report Style** option, *press  to toggle this option from Short to Long.* Then reset the monitor as described in **2.8 Settings Menu**.




Figure 14: Changing Report Style



Please use ESC/POS compatible serial interface thermal receipt printer and use null modem serial cable for connecting the thermal printer to the BP+. For more information, please contact Uscom.

### 2.13 Reviewing saved measurements

Measurements saved to the SD Card (if used) can be accessed through the Memory Mode menu for review.

Press . The Memory Mode menu will be displayed. If no SD Card is inserted, then the device will prompt you to insert one.

Use and to highlight the measurement to review. New measurements will be added to the end of the list. The screen can accommodate 40 measurements. The index number of the most recent measurement is displayed in a different colour. If less than 40 measurements are present on the SD Card, then '-' will be displayed in some positions.

Memory Mode			
1	2	3	4
5	6	7	8
9	10	11	12
13	14	15	16
17	18	19	20
21	22	23	24
25	26	27	28
29	30	31	32
33	34	35	36
37*	-	-	-
Selected 2		*Last Saved	

Figure 15: Selecting Measurement in Memory Mode

When more than 40 measurements are present on the SD Card, then the 40 measurements with the highest measurement numbers will be displayed by default when the SD Card is inserted. Otherwise, the Memory Mode will recall and display the last viewed page and cursor location whenever the Memory Mode menu is re-opened.

To view a lower measurement number that is not displayed on the current Memory Mode screen, use to move the cursor to select the first measurement index number on the screen and then press again to go to the previous page.

Vice versa, to view higher measurement numbers not displayed on the current Memory Mode screen, use to move the cursor to select the last measurement index number on the screen and then press again to go to the next page.



Memory Mode			
1	2	3	4
5	6	7	8
9	10	11	12
13	14	15	16
17	18	19	20
21	22	23	24
25	26	27	28
29	30	31	32
33	34	35	36
37	38	39	40
Press Start to select 40		*Last Saved	



Memory Mode			
41	42	43	44
45	46	47	48
49	50	51	52
53	54	55	56
57	58	59	60
61	62	63	64
65	66	67	68
69	70	71	72
73	74	75	76
77	78	79	80*

Press Start to select 41 \*Last Saved

Figure 16: Memory Mode Measurements Navigation

Press  to confirm your selection. The BP+ will begin to reprocess the measurement after which it will be displayed on the main screen. If you do not wish to confirm your selection, press .




Figure 17: Calculating Measurement in Memory Mode





Figure 18: BP Result Screen

## 3 Theory of Operation

### 3.1 Physiology of wave reflection



**Figure 19: Arteries of the upper body**

The Uscom BP+ monitor has been developed using Pulse Wave Oscillometric technology (referred to as **BP+**<sup>®</sup>) and a scientific understanding of pressure wave propagation in the arterial system. The arteries most relevant to BP+ measurement technique are shown in Figure 19.

The theory of wave reflections implies that the pressure at any location in the arterial tree can be considered to be the sum of forward- and backward-going pressure waves.


The initial forward-going pressure wave is generated by the contraction of the heart's left ventricle. Backward-going pressure waves are created when the forward-going pressure wave encounters a change in the properties and geometry within the arterial system. At such a point, the forward-going wave is partly reflected, creating the backward-going pressure wave. There are multiple reflection sites within the human arterial system.

### 3.2 Blood Pressure Measurement

The Cardiovascular monitor is designed to measure pressure-related information from the upper arm using an inflatable cuff. As shown in Figure 20, the pressure wave generated by the heart can be considered to take two paths to the upper arm cuff, generating what are known as the incident wave and reflected wave.


- The incident wave travels from the heart, through the ascending aorta, subclavian artery and brachial artery to reach the cuff.
- The reflected wave travels from the heart, through the ascending aorta, aortic arch and down the descending aorta to the effective reflection site in the abdominal aorta. At this location, some of the forward-going wave is reflected and travels back up the descending aorta, through the subclavian and brachial arteries before arriving at the cuff.

Uscom recommends measuring from the left arm. In cases where left-arm measurement is contraindicated, measurement from the right arm is possible. However, due to the additional complexity of the wave path to the right brachial artery, results from the right arm may not be comparable with left-arm results.



**Figure 20: Pressure wave reflection in the arterial system**

### 3.3 Display of Pulse Waveform and Rhythm strip



**Figure 21: BP Result Screen, mean beat feature points**

The following are the points of interest on the pulse waveform, as indicated in Figure 21.

- Start of the pulse,  $t_0, p_0$
- Peak of the incident wave,  $t_1, p_1$
- Trough between incident and reflected wave,  $t_2, p_2$
- Peak of the reflected wave,  $t_3, p_3$
- Trough of the dicrotic notch,  $t_4, p_4$
- Peak following the dicrotic notch,  $t_5, p_5$

### 3.4 Central Blood Pressure

The BP+ uses a physics-based model of the subclavian to brachial artery branch to calculate central blood pressure. This model has been validated against invasive pressure data obtained from the aortic arch, with measurements taken



on the left arm. Validation was performed in subjects aged between 30 and 80 years.

Central blood pressure measurement results may not be accurate in subjects who are not similar to those used in validation.

### 3.5 Pulse Wave Parameters Calculations

Referring to Figure 21, the equations for calculating pulse wave parameters displayed on the main screen and Pulse Wave Analysis screen are:

**Augmentation index (%):**

$$AI = \frac{p_3 - p_0}{p_1 - p_0}$$

**Pulse rate variability (milliseconds):**

$$PRV = \sqrt{\frac{1}{N-2} \sum_{n=2}^{N-1} (t_{0,n+1} - t_{0,n-1})^2}$$

**Pulse rate (beats per minute):**

$$\text{Pulse Wave PR} = \frac{1}{N-1} \sum_{n=2}^N (t_{0,n} - t_{0,n-1})$$

**Signal to noise ratio (decibels):**

$$SNR = 10 \log_{10} \left( \frac{\overline{(p(t) - \overline{p(t)})^2}}{\frac{1}{N} \sum_{n=1}^N (p_n(t) - p(t))^2} \right)$$

$$\text{where } \overline{f(t)} = \frac{1}{T} \int_{t=0}^T f(t) dt$$

Time to reflection (seconds):  $TR = t_3 - t_1$

Systolic ejection period (seconds):  $SEP = t_4 - t_0$

Reflected wave transit time (seconds):  $RWTT = t_2 - t_0$

Pulse wave pulse pressure (mmHg):  $PP = \max(p_1 - p_0, p_3 - p_0)$

Pulse wave pulse pressure variation (%):

$$DPP = \frac{\max_{n=1..N} (\max(p_{1,n} - p_{0,n}, p_{3,n} - p_{0,n})) - \min_{n=1..N} (\max(p_{1,n} - p_{0,n}, p_{3,n} - p_{0,n}))}{\frac{1}{N} \sum_{n=1}^N (\max(p_{1,n} - p_{0,n}, p_{3,n} - p_{0,n}))}$$

Maximum pulse wave pressure gradient (mmHg/s):

$$dP/dt(\max) = \max_t \left( \frac{dP(t)}{dt} \right)$$



## 4 Maintenance and Troubleshooting

### 4.1 Servicing



The Uscom BP+ contains no user-serviceable parts. Opening the BP+ enclosure voids any applicable warranty.

Some frequently asked questions are answered below. If your question is not answered, or for additional service information, please contact Uscom using the contact information at the front of this manual.

When contacting Uscom with service related enquiries, please provide the following information:

- Device Part Number
- Serial Number
- Description of any fault, including:
  - Circumstances and operation preceding the fault
  - Environment in which the device was being used at the time of the fault

### 4.2 Routine Maintenance

The BP+ is designed so that no periodic adjustment or calibration is required, if it is stored and used within the temperature and humidity ranges given in section 5 Specifications. The BP+ might not meet its performance specifications if stored or used outside of these ranges.

The accuracy of the cuff pressure measurement should be checked annually by a trained service technician.

The cuffs provided with the BP+ may become worn with routine use and require replacing.

### 4.3 Cleaning

The device may be cleaned using lint-free cloth dampened with a mild solution of detergent and water, 10% bleach solution, or a commercial disinfectant.

After cleaning, dry all areas with a lint-free cloth.



**Avoid ingress of liquid into any part of the BP+ or power supply unit. Do not submerge any component in liquid. This may cause fire or electrical shock.**

To clean the cuff and tubing, apply a mild detergent to a slightly damp cloth and wipe clean. Allow to dry thoroughly.

### 4.4 Disposal

Dispose or recycle the BP+ in accordance with regulations of your country.

Do not dispose of the BP+ as unsorted municipal waste at the end of the product's lifetime. The BP+ must be recycled in accordance with the WEEE Directive 2002/96/EC (for EU) or according to the regulations of your country. To arrange for return or disposal of the BP+ please contact your local supplier.

### 4.5 Frequently Asked Questions

**I cannot connect the air hose to the device.**

**I cannot connect the cuff to the extension hose.**

Push on the thumb button until you hear a click. Then try reinserting the coupling.

Sometimes the latching mechanism of the connector that usually locks the inserted coupling is activated without the coupling inserted. In this case it must be reset by pushing the thumb button.

**The cuff deflates during a measurement and “Safety Rule Violated” is displayed.**

There may have been an overpressure condition. For safety reasons, the cuff will deflate if the device detects a very high pressure in the cuff, or if a significant pressure is held for a long period of time. Turn the device off and then on again to reset the safety lock.

**The screen is completely white when I turn the device on.**

**The screen flickers when I turn the device on.**

**The device resets by itself.**

Check that an appropriate AC adapter is being used.

**I get the message “Cuff Too Tight or Kinked Hose” when trying to take a measurement.**

The device is having trouble inflating the cuff. Ensure that the cuff is on firmly but not tightly, and that the hose to the cuff is not kinked, compressed or blocked.

**I get the message “Invalid Signal” after taking a measurement.**

This can occur when the signal is too noisy to process. Ensure the cuff is applied firmly and the subject remains still during the measurement.

**I get the message “Cuff Not Detected” when trying to take a measurement.**

The cuff is not properly attached. Check all pneumatic connections are securely latched. In rare cases the O-ring on the connector may be damaged or missing. Contact Uscom for service information.

**The signal is consistently very noisy (SNR < 3) even under ideal measurement conditions.**

There may be a slow leak in the pneumatic circuit. Check for damaged air connectors (particularly the O-ring on each connector), split hoses or punctured cuffs. You may wish to use a different hose and cuff set to verify the problem.

If the problem persists, contact Uscom.

**The measured blood pressure is unexpected.**

**The measurements change a lot each time they are measured.**

Ensure that good blood pressure measurement practice is being followed as described in Section 2.1 Preparing for a measurement. Multiple successive measurements may be used to calculate an average value for all parameters.

## 4.6 Error Message

### 4.6.1 Main Screen Message

Message	Reason	User Action
<b>Connection Failed Or No Connection</b>	Internal fault condition.	Consult Uscom / Uscom certified servicing technician
	Wrong power adapter used.	Check the power adapter being used is the one provided with the device, as specified in Section <b>5 Specification</b>
<b>Measurement Timeout</b>	Measurement could not be completed within the maximum permitted time. Safety system has cancelled the measurement.	Redo another measurement. Ask the patient to relax, sit still, not to talk, no arm or finger movement on the arm being measured (left), etc. Refer to section <b>2.1 &amp; 2.2</b> .
<b>Systolic BP too High</b>	The upper arm systolic pressure measured is too high. The device cannot proceed to do the Pulse Wave measurement	If error persists, consult Uscom / certified servicing technician (Please make a note of any error number displayed along with the error message.)
<b>Unable to Calculate Results Unable to Find Features Unable to Find Pulse Wave</b>	Calculation error	
<b>Unable to Reinflate</b>	The inflation target for doing the Pulse Wave measurement is too high	
<b>Unable to Start</b>	Device detects an internal error when it starts the blood pressure measurement	
<b>Invalid Signal</b>	The signal noise ratio is too low and out of the valid range.	
<b>Poor Signal</b>	The signal quality is poor	
<b>Quality too low to Auto Save</b>	The signal quality is too low and the measurement will not be saved	
<b>Unable to measure BP: Please Repeat (####)</b> (Where #### is the error number)	Device detects a blood pressure measurement related error when it is doing the blood pressure measurement	
<b>Unable to measure BP: Timeout (####)</b> (Where #### is the error number)	Blood pressure measurement Timeout	
<b>Unable to measure BP: (####)</b> (Where #### is the error number)	Blood pressure measurement failed	Consult Uscom/Uscom certified servicing technician with the BP error number

Message	Reason	User Action
<b>Unable to measure BP: Check Pneumatics (####)</b> (Where #### is the error number)	Device detects pneumatics error when it is doing the blood pressure measurement	Check the cuff and hose is attached correctly. If the problem persist, consult Uscom or Uscom certified servicing technician. Please make a note of the error number.
<b>Invalid Measurement Data</b>	The target measurement record being retrieved from the SD card is corrupted	Delete the corrupted record from the SD card. This can be done by inserting the SD card into a computer and deleting the appropriate file.
<b>Measurement not saved</b>	No SD inserted when the device is trying to save the measurement data.	Insert a compatible SD card
	When the device is trying to save the measurement data and an oversize SD card is detected. The device does not support a SD card with a capacity larger than 2GB	Check the capacity of the SD card and make sure it is not larger than 2GB
	When device is trying to save the measurement data and the SD card is full	Delete older record(s) or replace with an empty SD card.
	When device is trying to save the measurement data and the SD card is not formatted correctly.	Format the SD card to FAT16 or FAT32 on a PC.
<b>Overpressure</b>	Pressure is over the limit	Remove the cuff from the arm and device. <u>Note:</u> Device will automatically deflate the cuff pressure immediately.
<b>Safety Rule Violated</b>		

#### 4.6.2 Memory Mode Message

Message	Reason	User Action
<b>Please replace SD Card. Full or Invalid</b>	SD card is full	Delete the old record(s) or replace with an empty SD card
<b>Unsupported SD Card</b>	Invalid SD card format	Format the SD card to FAT16 or FAT32 on computer
	Large capacity of SD card is inserted, device does not support capacity larger than 2GB	Check the capacity of the SD card and make sure it is not larger than 2GB
<b>Slot is Empty</b>	No SD card is inserted	Insert the compatible SD card

### 4.6.3 SD Card Icon Message

Message	Reason	User Action
<b>Not Saved</b>	No SD inserted when device is trying to save the measurement data.	Insert a compatible SD card
	when device is trying to save the measurement data: - Oversize SD card is detected, device does not support capacity larger than 2GB	Check the capacity of the SD card and make sure it is not larger than 2GB
	When device is trying to save the measurement data. - SD card is full	Delete the old record(s) or replace with an empty SD card
	When device is trying to save the measurement data. - The format of SD card is detected to be invalid	Format the SD card to FAT16 or FAT32
<b>Save Failed</b>	The measurement will not be saved due to the calculation error	Redo another measurement and make sure sit still and no arm movement, refer to section <b>2.1 and 2.2</b> . If error persists, consult Uscom / certified servicing technician.

## 5 Specifications

### Non-invasive upper-arm blood pressure

- Oscillometric method using step controlled deflation
- Adult and pediatric patients
- BP Measurement range:
  - Systolic: 60 to 245 mmHg
  - Diastolic: 30 to 160 mmHg
  - Mean: 40 to 190 mmHg
  - Cuff Pressure: 0 to 300 mmHg
- BP Measurement accuracy:
  - $\leq \pm 3$  mmHg over measurement range
- Pulse Rate range: 35 to 199 bpm
- Pulse Rate accuracy:  $\pm 5\%$
- Power-on self-test
- Typical measurement time within 60 seconds

### Non-invasive central blood pressure

- Model-based transform of peripheral waveform reports Systolic and Diastolic central pressure

### Pulse wave measurement

- Arterial measures:
  - Peripheral Augmentation Index (AI)
  - Time of Reflection (TR)
  - Reflected Wave Transit Time (RWTT)
- Cardiac function measures:
  - Oscillometric pulse pressure (PP)
  - Maximum pulse wave pressure gradient ( $dP/dt|_{max}$ )
  - Pulse wave pulse pressure variation (DPP)
  - Systolic ejection period (SEP)
  - Pulse rate variability (PRV)
- Signal to Noise ratio (SNR)
- Data Sampling Rate: 200Hz

### Physical dimensions

- Size: 156 mm x 157 mm x 119 mm
- Weight: 500 g

### Display and memory

- 320x240 pixel matrix 3.5" TFT colour display
- Pulse wave waveform display
- Graphical rhythm strip display
- Status and error display
- SD Card slot, hot swappable with cards up to 2 GB

### Blood pressure cuff and hose

- Lifetime antimicrobial treatment
- Small Adult, Adult, Extra Long Adult and Large Adult cuffs available
- 1.2 m extension hose

### Computer interface

- Computer interface specification (serial port) available for customer configuration

### Power supply

- Power Supply Unit: GlobTek Inc model GTM41060-1706 medical grade, AC/DC wall adapter
- AC Input: 100–240 V<sub>AC</sub>, 50–60 Hz, 0.6 A
- DC Output: 6 V<sub>DC</sub>, 2.8 A
- Use of the BP+ with the Power Supply Unit is considered an ME system

### Environmental specifications

- Operating Temperature: 10 to 40 °C (50 to 104 °F)
- Operating Relative Humidity: 15 to 80% (non-condensing)
- Operating Atmospheric Pressure: 80 to 103 kPa
- Transport and Storage Temperature: -20 to 50 °C (-4 to 122 °F)
- Transport and Storage Relative Humidity: 15 to 80% (non-condensing)

### IEC 60601-1 standard classification

- Class II Equipment
- Type BF Applied Part (the blood pressure cuff is considered the applied part)
- Continuously powered with measurements being taken upon demand
- Equipment not suitable for use in the presence of a flammable anesthetic mixture, with air or oxygen or nitrous oxide

### Warranty

- 12 month warranty on the device
- 6 month warranty on cuffs and extension hose

## 5.1 Electromagnetic Compatibility

Guidance and manufacturer's declaration - electromagnetic emissions			
The BP+ is intended for use in the electromagnetic environment specified below. The customer or the user of the BP+ should assure that it is used in such an environment.			
Emissions test	Compliance	Electromagnetic environment - guidance	
RF emissions CISPR 11	Group 1	The BP+ uses RF energy only for its internal function. Therefore, its RF emissions are very low and not likely to cause any interference in nearby electronic equipment. The BP+ is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.	
RF emissions CISPR 11	Class B		
Harmonic emissions IEC 61000-3-2	Class A		
Voltage fluctuations / flicker emissions IEC 61000-3-3	Complies		
The BP+ should not be used adjacent to or stacked with other equipment. If adjacent or stacked use is necessary, the BP+ should be observed to verify normal operation in the configuration in which it will be used.			
Guidance and manufacturer's declaration - electromagnetic immunity			
The BP+ is intended for use in the electromagnetic environment specified below. The customer or the user of the BP+ should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge, ESD (IEC 61000-4-2)	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst (IEC 61000-4-4)	±2 kV for power supply lines ±1 kV for input / output lines	±2 kV for power supply lines ±1 kV for input / output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge (IEC 61000-4-5)	±1 kV line(s) to line(s) ±2 kV line(s) to earth	±1 kV line(s) to line(s) ±2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines (IEC 61000-4-11)	<5 % $U_T$ (>95 % dip in $U_T$ ) for 0.5 cycle 40 % $U_T$ (60% dip in $U_T$ ) for 5 cycles 70 % $U_T$ (30 % dip in $U_T$ ) for 25 cycles <5 % $U_T$ (>95 % dip in $U_T$ ) for 5 sec	<5 % $U_T$ (>95 % dip in $U_T$ ) for 0.5 cycle 40 % $U_T$ (60% dip in $U_T$ ) for 5 cycles 70 % $U_T$ (30 % dip in $U_T$ ) for 25 cycles <5 % $U_T$ (>95 % dip in $U_T$ ) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the BP+ requires continued operation during power mains interruptions, it is recommended that the BP+ be powered from an uninterruptible power supply.
Power frequency (50/60 Hz) magnetic field (IEC 61000-4-8)	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE $U_T$ is the AC mains voltage prior to application of the test level.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF (IEC 61000-4-6)	3 Vrms 150 kHz to 80 MHz	3 V	Portable and mobile RF communications equipment should be used no closer to any part of the BP+, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = 1.2 \sqrt{P}$ $d = 1.2 \sqrt{P}$ 80 MHz to 800 MHz $d = 2.3 \sqrt{P}$ 800 MHz to 2.5 GHz where $P$ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and $d$ is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, <sup>a</sup> should be less than the compliance level in each frequency range. <sup>b</sup> Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF (IEC 61000-4-3)	3 V/m 80 MHz to 2.5 GHz	3 V/m	
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
<sup>a</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular / cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the BP+ is used exceeds the applicable RF compliance level above, the BP+ should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the BP+.			
<sup>b</sup> Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			
Recommended separation distances between portable and mobile RF communications equipment and the BP+			
The BP+ is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the BP+ can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the BP+ as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 1.2 \sqrt{P}$	80 MHz to 800 MHz $d = 1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance $d$ in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where $P$ is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
<b>Essential Performance of the BP+ as required by IEC 60601-1-2 § 3.201.2 is defined as 1) opening of cuff pressure valves if detected cuff pressure exceeds 300 mmHg or 15 mmHg for more than 3 minutes; 2) no degradation of blood pressure measurement accuracy beyond that given in the specification.</b>			

## 6 Accessories and Spare Parts

### Adult Extra Long Cuff

Welch Allyn Part: REUSE-11L

Part: U50304



This is a re-useable, latex free cuff suitable for arm circumference 25 cm to 34 cm. It will also fit arms of greater circumference. This is the most suitable cuff for general adult use.

### Small Adult Cuff

Welch Allyn Part: REUSE-10

Part: U50303



This is a re-useable, latex free cuff suitable for arm circumference 20 to 26 cm.

### Large Adult Cuff

Welch Allyn Part: REUSE-12

Part: U50305



This is a re-useable, latex free cuff suitable for arm circumference 32 cm to 43cm.

### Extension Hose

Part: U50301



This is a 1.2 m latex free extension hose to connect between the device and the cuff.

### Medical grade mains power adapter with international blade set

Model Number: GlobTek Inc GTM41060-1706

Part: U50307



Suitable for use as supply of external DC power to the BP+. Comes with US, Europe, United Kingdom, and Australia/New Zealand blade set. Has medical grade certification for use in hospital environments.